UPUTSTVO ZA UPOTREBU

(SRB)

XLD - modifikovani Agar Plate

Selektivna podloga za izolaciju i enumeraciju *Salmonella Typhi* i drugih *Salmonella* vrsta.

Sadržaj pakovanja:

Šifra artikla (pakovanja)	Opis	Šifra primarnog pakovanja:	Broj podloga
PRM031IV20	Podloga izlivena u petri posudama od ø90	PRM031I	20
PRM031IV60			60
PRM031IV240			240
PRM031IM40	Podloga izlivena u petri posudama od ø50		40

Uputstva

Pod aseptičnim uslovima inokulisati ploču metodom površinskog zasejavanja. Nakon inkubacije posmatrati rast i boju kolonija.

Princip i interpretacija

Xylose Lysine Deoxycholate (XLD) agar je napravio Taylor (1-5) za izolaciju i diferencijaciju enteričnih patogena uključujući i *Salmonella Typhi* od drugih vrsta roda *Salmonella*. XLD Agar, modifikovani se preporučuje za selektivnu izolaciju i enumeraciju *Salmonella Typhi* i drugih *Salmonella* vrsta u skladu sa ISO Committee (20).

XLD agar je podloga pojačane selektivnosti i osetljivosti u poređenju sa drugim podlogama kao što su: SS Agar, EMB Agar i Bizmut-sulfitni agar (2,4,6 i 13-16). Formulacija podloge onemogućava prerastanje *Salmonella* spp. i *Shigella* spp. od strane drugih mikroorganizama (17).

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Literatura

- Taylor W. L., 1965, Am. J. Clin. Pathol., 44:471-475.
- Taylor W. L. and Harris B., 1965, Am. J. Clin. Pathol., 44:476.
- Taylor W. L. and Harris B., 1967, Am. J. Clin. Pathol., 48:350.
- Taylor W. L. and Schelhart B., 1967, Am. J. Clin. Pathol., 48:356.
- Taylor W. L. and Schelhart B., 1968, Am. J. Clin. Pathol., 16:1387.
- Taylor W. L. and Schelhart B., 1969, Appl. Microbiol., 18:393-395.
- Chadwick P., Delisle G. H and Byer M., 1974, Can. J. Microbiol., 20, 1653-1664.
- Downes F. P. and Ito K. (Ed.), 2001, Compendium of Methods for the Microbiological examination of Foods, 4th Ed., APHA Inc. Washington D.C.
- Eaton A. D., Clesceri L. S., Rice E. W., and Greenberg A. W., (Eds.), 2005, Standard Methods for the Examination of Water and Wastewater, 21st Ed., APHA, Washington, D.C.
- Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
- Williams H., (Ed.), 2005, Official Methods of Analysis of the Association of Official Analytical Chemists, 19th Ed., AOAC, Washington, D.C.
- FDA Bacteriological Analytical Manual, 2005, 18th Ed., AOAC, Washington, D.C.
- Dunn C. and Martin W. J., 1971, Appl. Microbiol., 22, 17-22.
- Rollender M. A., Beckford O., Belsky R. D and Kostroff B. 1969, Am. J. Clin. Pathol., 51, 284-286.
- Taylor W. L. and Schelhart B., 1969, Appl. Micro. 18, 1387-1392.
- MacCarthy M. D., 1966, N. Z. J. Med. Lab. Technol., 20, 127-131.
- Isenberg H. D., Kominos S., and Sigal M., 1969, Appl Microbiol., 18, 656-659.
- Aspinall S. T., Hindle M. A. and Hutchinson D. N., 1992, 19. Eur. J. Clin. Microbiol., Inf. Dis. 11, 936-939. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore.
- International Organization for Standardization (ISO), 2002, Draft ISO/DIS 6579:2002.

Proizvođač: ProMedia doo, Kralja Petra I 114, 23300 Kikinda, Srbija

Mesto proizvodnje: PC - PROIZVODNJA PRIPREMLJENIH MIKROBIOLOŠKIH HRANLJIVIH MEDIJUMA PROREADY, KIKINDA.

Put za pristanište 140, 23300 Kikinda, Srbija; Tel: +381 230/315-045; mail: office@proready.rs; www.proready.rs

INSTRUCTION FOR USE

(EN)

XLD Agar- Modified Agar Plate

Media for selective isolation and enumeration of *Salmonella* Typhi and other *Salmonella* species.

Package contents:

Item code (packaging) REF	Description	Primary packaging code:	Number of products
PRM031IV20	Substrate poured into petri dishes of ø90	PRM0311	20
PRM031IV60			60
PRM031IV240			240
PRM031IM40	Substrate poured into petri dishes of ø50		40

Directions

Surface spread the test inoculum aseptically on the plate. After incubation, observe growth and color of colonies.

Principle And Interpretation

XLD Agar was formulated by Taylor (4-8) for the isolation and differentiation of enteric pathogens including *Salmonella* Typhi from other *Salmonella* species. XLD Agar, Modified is recommended for selective isolation and enumeration of *Salmonella* Typhi from other *Salmonella* species in accordance with ISO Committee (3).

XLD agar has enhanced selectivity and sensitivity compared to other media such as: SS Agar, EMB Agar and Bismuth sulphite agar. Its formulation prevents the overgrowth of *Salmonella* spp. and *Shigella* spp. by other microorganisms (2).

This medium contains yeast extract, which provides nitrogen and vitamins required for growth. Xylose, lactose and sucrose are sources of fermentable carbohydrates, however xylose is mainly incorporated into the medium since it is not fermented by *Shigella* spp. but practically by all enterobacteria. This helps in the differentiation of *Shigella* species. Sodium chloride maintains the osmotic balance of the medium. Lysine is included to differentiate the *Salmonella* group from the non-pathogens. *Salmonella* spp rapidly ferment xylose and exhaust carbon source. Subsequently lysine is decarboxylated by the enzyme lysine decarboxylase to form amines with reversion to an alkaline pH that mimics the *Shigella* spp. reaction. However, to prevent this reaction by lysine-positive coliforms, lactose and sucrose are added to produce acid in excess. Degradation of xylose, lactose and sucrose to acid causes phenol red indicator to change its colour to yellow. Bacteria that decarboxylate lysine to cadaverine can be recognized by the appearance of a red colouration around the colonies due to an increase in pH. These reactions can proceed simultaneously or successively, and this may cause the pH indicator to exhibit various shades of colour or it may change its colour from yellow to red on prolonged incubation.

In order to ensure the differentiation of colonies, sodium thiosulphate and ferric ammonium citrate are included in the medium as H₂S-indicators, which provide visualization of H₂S production, which is manifested by the formation of colonies with a black center. Non-pathogenic bacteria that produce H₂S do not decarboxylate lysine, which is why the production of acids in their case prevents blackening of colonies (1).

XLD agar is both a differential and a selective medium. Sodium deoxycholate is a selective agent as an inhibitor of Gram-positive microorganisms.

Quality control

The data and results of quality control are given in the certificate of analysis for each lot.

Storage and shelf life

Storage between 15-25°C. After opening store between 2-8°C. Use before expiry date on the label.

Warning and precautions

This product does not contain hazardous substances in concentrations that are above the prescribed limits set by applicable legislation and are therefore is not classified as hazardous. However, it is recommended to follow the guidelines provided in the safety data sheet for proper use. This product is intended for laboratory use only by a professionally trained person.

Do not use the product if the primary packaging is damaged or the product does not meet the stated characteristics.

Disposal

Waste disposal must be in accordance with national and local regulations. Each laboratory is responsible for handling and disposing of waste generated during operation.

Symbols used on labels

	This side up		Catalogue number
	Do not expose directly to sunlight		Batch code
	Consult instructions for use		Use-by date
	Do not re-use		Temperature limit
	Pack size		Manufacturer

Reference

- Taylor W. L., 1965, Am. J. Clin. Pathol., 44:471-475.
- Taylor W. L. and Harris B., 1965, Am. J. Clin. Pathol., 44:476.
- Taylor W. L. and Harris B., 1967, Am. J. Clin. Pathol., 48:350.
- Taylor W. L. and Schelhart B., 1967, Am. J. Clin. Pathol., 48:356.
- Taylor W. L. and Schelhart B., 1968, Am. J. Clin. Pathol., 16:1387.
- Taylor W. L. and Schelhart B., 1969, Appl. Microbiol., 18:393-395.
- Chadwick P., Delisle G. H and Byer M., 1974, Can. J. Microbiol., 20, 1653-1664.
- Downes F. P. and Ito K. (Ed.), 2001, Compendium of Methods for the Microbiological examination of Foods, 4th Ed., APHA Inc. Washington D.C.
- Eaton A. D., Clesceri L. S., Rice E. W., and Greenberg A. W., (Eds.), 2005, Standard Methods for the Examination of Water and Wastewater, 21st Ed., APHA, Washington, D.C.
- Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
- Williams H., (Ed.), 2005, Official Methods of Analysis of the Association of Official Analytical Chemists, 19th Ed., AOAC, Washington, D.C.
- FDA Bacteriological Analytical Manual, 2005, 18th Ed., AOAC, Washington, D.C.
- Dunn C. and Martin W. J., 1971, Appl. Microbiol., 22, 17-22.
- Rollender M. A., Beckford O., Belsky R. D and Kostroff B. 1969, Am. J. Clin. Pathol., 51, 284-286.
- Taylor W. L. and Schelhart B., 1969, Appl. Micro. 18, 1387-1392.
- MacCarthy M. D., 1966, N. Z. J. Med. Lab. Technol., 20, 127-131.
- Isenberg H. D., Kominos S., and Sigel M., 1969, Appl Microbiol., 18, 656-659.
- Aspinall S. T., Hindle M. A. and Hutchinson D. N., 1992, 19. Eur. J. Clin. Microbiol., Inf. Dis. 11, 936-939. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore.
- International Organization for Standardization (ISO), 2002, Draft ISO/DIS 6579:2002.