UPUTSTVO ZA UPOTREBU (SRB)

Sabouraud Dextrose Agar Plate

Podloga za kultivaciju kvasaca, plesni i acidodurnih bakterija.

Sadržaj pakovanja:

Šifra artikla (pakovanja)	Opis	Šifra primarnog pakovanja:	Broj podloga
PRM063V20	Podloga izlivena u petri posudama od ø90	PRM063	20
PRM063V60			60
PRM063V240			240
PRM063M40	Podloga izlivena u petri posudama od ø50		40

Uputstva

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


(EN)

Sabouraud Maltose Agar Plate

Medium is used for the cultivation of yeasts, moulds and aciduric bacteria.

Package contents:

Item code (packaging) REF	Description	Primary packaging code:	Number of products
PRM063V20	Substrate poured into petri dishes of ø90	PRM063	20
PRM063V60			60
PRM063V240			240
PRM063M40	Substrate poured into petri dishes of ø50		40

Directions

Surface spread the test inoculum aseptically on the plate.

Principle and interpretation

Sabouraud Dextrose Agar is Carliers modification (1) of the formulation described by Sabouraud (2) for the cultivation of fungi (yeasts, moulds), particularly useful for the fungi associated with skin infections. This medium is also employed to determine microbial contamination in food, cosmetics and clinical specimens (3).

Mycological peptone provides nitrogenous compounds. Dextrose provides an energy source. High dextrose concentration and low pH favours fungal growth and inhibits contaminating bacteria from test samples (4).

Some pathogenic fungi may produce infective spores which are easily dispersed in air, so examination should be carried out in safety cabinet.

Quality control

The data and results of quality control are given in the certificate of analysis for each lot.

Storage and shelf life

Storage between 15-25°C. After opening storage between 2-8°C . Use before expiry date on the label.



Warning and precautions


In vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Symbols used on labels

	European Conformity mark		This side up
	is an in vitro diagnostic medical device (IVD)		Catalogue number
	Do not expose directly to sunlight		Batch code
	Consult instructions for use		Use-by date
	Do not re-use		Temperature limit
	Pack size		Manufacturer
	European Authorized Representative (Authorised Representative)		

	Salus Cons kft. 6722 Szeged, Bécsi krt 23, HUNGARY e-mail: office@saluscons.com
---	--

Reference

1. Carlier G. I. M., 1948, Brit. J. Derm. Syph., 60:61.
2. Sabouraud K., 1892, Ann. Dermatol. Syphilol, 3:1061.
3. Bacteriological Analytical Manual, 8th Edition, Revision A, 1998. AOAC, Washington D.C.
4. Murray PR, Baren EJ, Jorgensen JH, Pfaller MA, Tenover FC, Tenover FC (editors) 2003, Manual of clinical Microbiology, 8th ed., ASM, Washington, D.C.