UPUTSTVO ZA UPOTREBU

(SRB)

Sheep Blood Agar Plate

Podloga za izolaciju i kultivaciju visoko zahtevnih patogenih mikroorganizama kao što su Neisseria, streptokoke itd.

Sadržaj pakovanja:

Šifra artikla (pakovanja) REF	Opis	Šifra primarnog pakovanja:	Broj podloga
PRM073V20	Podloga izlivena u petri posudama od ø90	PRM073	20
PRM073V60			60
PRM073V240			240
PRM073M40	Podloga izlivena u petri posudama od ø50		40

Uputstva

Pod aseptičnim uslovima inokulisati ploču metodom površinskog zasejavanja. Nakon inkubacije posmatrati rast i hemolizu oko kolonija.

Princip i interpretacija

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(EN)

Sheep Blood Agar Plate

Medium is used for isolation and cultivation of fastidious pathogenic microorganisms like Neisseria, Streptococci etc.

Package contents:

Item code (packaging) REF	Description	Primary packaging code:	Number of products
PRM073V20	Substrate poured into petri dishes of ø90	PRM073	20
PRM073V60			60
PRM073V240			240
PRM073M40	Substrate poured into petri dishes of ø50		40

Directions

Surface spread the test inoculum aseptically on the plate. After incubation, observe growth and color of colonies.

Principle and interpretation

Sheep Blood Agar is a highly nutritious medium used for isolation and cultivation of fastidious microorganisms. Sheep Blood Agar is recommended by APHA (5) and Standard Methods (6, 7) for testing of food samples. Beef extract and tryptose provides carbon, nitrogen, amino acids and vitamins. Sodium chloride helps in maintaining the osmotic equilibrium of the medium. Presence of blood makes the medium more nutritious by providing additional growth factors required by fastidious organisms. It also helps in visualizing the haemolytic reactions. However, haemolytic reactions depend on the animal blood used. Sheep blood gives best results for Group A Streptococci (8), but sheep blood fails to support growth of Haemophilus haemolyticus since sheep blood is deficient in pyridine nucleotides..

Quality control

The data and results of quality control are given in the certificate of analysis for each lot.

Storage and shelf life

Storage between 2-8°C. Use before expiry date on the label.

Warning and precautions

In vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Symbols used on labels

CE	European Conformity mark	↑↑	This side up
IVD	is an in vitro diagnostic medical device (IVD)	REF	Catalogue number
	Do not expose directly to sunlight	LOT	Batch code
	Consult instructions for use		Use-by date
	Do not re-use		Temperature limit
	Pack size		Manufacturer
EC REP	European Authorized Representative (Authorised Representative)		

EC REP	Salus Cons kft. 6722 Szeged, Bécsi krt 23, HUNGARY e-mail: office@saluscons.com
--------	--

Reference

1. Noble W. C., 1962, J. Clin. Pathol., 15:552.
2. Hansen N. H., 1962, J. Appl. Bacteriol., 25:46.
3. Hayes P. R., 1963, J. Gen. Microbiol., 30:1.
4. Schuber J. H., Edwards P. R. and Ramsere C. H., 1959, J. Bacteriol., 77:648.
5. Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, DC.
6. U.S. Food and Drug Administration, 1995, Bacteriological Analytical Manual, 8th Ed., AOAC International, Gaithersburg, Md.
7. Atlas R. M., 1993, Handbook of Microbiology of Microbiological Media, CRC Press, Boca Raton, Fla.
8. Snavey J. G. and Brahier J., 1960, Am. J. Clin. Pathol., 33:511.
9. Murray P. R., Baron J. H., Jorgensen J. H., Pfaller M. A., Tenover J. C., (Eds.), 8th Ed., 2003, Manual of Clinical Microbiology, ASM, Washington, D.C.