UPUTSTVO ZA UPOTREBU (SRB)

Mueller Hinton Agar Plate w/ Horse Blood

Podloga za testiranje antimikrobne osetljivosti zahtevnih organizama (pneumococcus i drugih streptococci, Haemophilus spp. i Moraxella spp.) u skladu sa smernicama EUCAST-a.

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


Literatura

- Murray P. R., Baron J. H., Pfaller M. A., Tenover J. C. and Tenover F. C., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
- Mueller J. H. and Hinton J., 1941, Proc. Soc. Exp. Biol. Med., 48:330.
- National Committee for Clinical Laboratory Standards, 2000, Approved Standard: M7-A5. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that grow aerobically, 5th Ed., NCCLS, Wayne, Pa.
- NCCLS Approved Standard: ASM-2, 1979, Performance Standards for Antimicrobial disc Susceptibility Tests, 2nd Ed., National Committee for Clin. Lab. Standards.
- Bauer A. W., Kirby W. M., Sherris J. L. and Turck M., 1966, Am. J. Clin. Pathol., 45:493.
- Present Status and Future Work, WHO Sponsored collaborative study, Chicago, Oct. 1967.
- Ericsson H. M. and Sherris J. L., 1971, Acta Pathol. Microbiol., Scand. Sect B Suppl., 217:1.
- National Committee for Clinical Laboratory Standards, 1986, Proposed Standards, M6-P, NCCLS, Villanova, Pa.
- Koch A. E. and Burchall J. J., 1971, Appl. Microbiol., 22: 812.
- Ferone R. Bushby R. M., Burchall J. J., Moore W. D., Smith D., 1975, Antimicrob. Agents chemotherap., 7 : 91.
- Pollock H. M., Minshew B. H., Kenney M. A., Schoenkecht F. D., 1978 , Antimicrob. Agents Chemotherap.; 14:360.
- DAmato R. F., and Thornsberry C., 1979, Curr. Microbiol., 2 : 135.

Broj rešenja o registraciji: 515-02-02534-22-003

Proizvođač: ProMedia doo, Kralja Petra I 114, 23300 Kikinda, Srbija

Mesto proizvodnje: PC - PROIZVODNJA PRIPREMLJENIH MIKROBIOLOŠKIH HRANLJVIH MEDIJUMA PROREADY, KIKINDA.

Put za pristanište 140, 23300 Kikinda, Srbija; Tel: +381 230/315-045; mail: office@proready.rs; www.proready.rs