UPUTSTVO ZA UPOTREBU (SRB)

Chromogenic UTI Agar Plate

Podloga za identifikaciju, diferencijaciju i potvrdu enteričnih bakterija iz uzoraka kao što je urin, koji mogu da sadrže veliki broj *Proteus* vrsta, kao i potencijalno patogene Gram-pozitivne organizme.

Sadržaj pakovanja:

Šifra artikla (pakovanja) REF	Opis	Šifra primarnog pakovanja:	Broj podloga
PRM1418V20	Podloga izlivena u petri posudama od ø90	PRM1418	20
PRM1418V60			60
PRM1418V240			240
PRM1418M40	Podloga izlivena u petri posudama od ø50		40

Uputstva

Pod aseptičnim uslovima inokulisati ploču metodom površinskog zasejavanja. Nakon inkubacije posmatrati rast i boju kolonija.

Princip i interpretacija

Chromogenic UTI Agar je formulisan na osnovu rada koji su sproveli Pezzlo (1), Wilkie i sar. (2), Friedman i sar. (3), Murray i sar. (4), Soriano i Ponte (5) i Merlino i sar. (6). Ova podloga se može koristiti umesto MacConkey Agara za izolaciju i potvrdu različitih organizama. Omogućava brzu identifikaciju nekih Gram-negativnih bakterija i nekih Gram-pozitivnih bakterija na osnovu različitih boja kolonija proizvedenih reakcijom enzima specifičnih za rodove ili vrste sa dva hromogena supstrata.

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


(EN)

Chromogenic UTI Agar Plate

Medium for identification, differentiation and confirmation of enteric bacteria from specimens such as urine which may contain large number of Proteus species as well as potentially pathogenic Gram-positive organisms.

Package contents:

Item code (packaging) REF	Description	Primary packaging code:	Number of products
PRM1418V20	Substrate poured into petri dishes of \varnothing 90	PRM1418	20
PRM1418V60			60
PRM1418V240			240
PRM1418M40	Substrate poured into petri dishes of \varnothing 50		40

Directions

Surface spread the test inoculum aseptically on the plate. After incubation, observe growth and color of colonies.

Principle and interpretation

Chromogenic UTI Agar is formulated on the basis of work carried out by Pezzlo (1), Wilkie et al (2), Friedman et al (3), Murray et al (4), Soriano and Ponte (5) and Merlino et al (6). This media is used instead of MacConkey Agar for isolation and confirmation of various microorganisms. It facilitates and expedites the identification of some Gram-negative bacteria and some Gram-positive bacteria on the basis of different contrasted colony colours produced by reactions of genus or species specific enzymes with two chromogenic substrates.

Enzymes produced by Enterococcus species, Escherichia coli and coliforms cleave the chromogenic substrates incorporated in the medium. Presence of rich source of phenylalanine and tryptophan from peptone and casein enzymic hydrolysate provides an indication of tryptophan deaminase activity, revealed with TDA Reagent (R036) indicating the presence of Proteus species, Morganella species and Providencia species, which appear brown. One chromogenic substrate is cleaved by β -glucosidase possessed by Enterococci resulting in formation of blue colonies. E. coli produce purple-magenta colonies due to the enzyme β -D-galactosidase which cleaves the other chromogenic substrate. Further confirmation of E.coli can be done by performing indole test using DMACA Reagent (R035). Also, some strains of Enterobacter cloacae lacking β -glucosidase show pink-colonies indistinguishable from E.coli. The DMACA reagent for indole test (should be performed on filter paper) distinguishes between E.coli and Enterobacter, and also between Proteus mirabilis and other species. Coliforms produce purple coloured colonies due to cleavage of both the chromogenic substrates.

Peptone, casein enzymic hydrolysate and meat extract B# provides nitrogenous, carbonaceous compounds and other essential growth nutrients.

Quality control

The data and results of quality control are given in the certificate of analysis for each lot.

Storage and shelf life

Storage between 2-8°C. Use before expiry date on the label.


Warning and precautions

In vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Symbols used on labels

	European Conformity mark		This side up
	is an in vitro diagnostic medical device (IVD)		Catalogue number
	Do not expose directly to sunlight		Batch code
	Consult instructions for use		Use-by date
	Do not re-use		Temperature limit
	Pack size		Manufacturer
	European Authorized Representative (Authorised Representative)		

	Salus Cons kft. 6722 Szeged, Bécsi krt 23, HUNGARY e-mail: office@saluscons.com
---	--

Reference

1. Pezzlo M., (1998), Clinical Microbiology Reviews, 1:268-280.
2. Wilkie M. E., Almond M. K. and Marsh F. p., (1992), British Medical Journal, 305:1137-1141.
3. Friedman M. P. et al. (1991), Journal of Clinical Microbiology, 29:2385-2389.
4. Murray P., Traynor P. And Hopson D., (1992), Journal of Clinical Microbiology, 30:1600-1601.
5. Soriano F. and Ponte C., (1992), Journal of Clinical Microbiology, 30:3033-3034.
6. Merlino et al. (1995), Abstr. Austr. Microbiol., 16(4):17-3.