UPUTSTVO ZA UPOTREBU (SRB)

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
Literatura

- Anthony BF, Okada DM, Hobel CJ. Epidemiology of group B Streptococcus: longitudinal observations during pregnancy. J Infect Dis 1978; 137:524-30.
- Murray P.R., Baron J.H., Tenover FC, Tenover FC, (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
- Prevention of perinatal group B Streptococcal disease: a public health perspective. Centres for Disease control and Prevention. MMWR Recomm Rep 1996; 51:1-22
- NHS Processing swabs for Group B Streptococcal carriage Issue no.2.1,2006

Broj rešenja o registraciji: 515-02-02534-22-003

INSTRUCTION FOR USE (EN)

Chromogenic Strep B Selective Agar Plate

Medium is recommended for selective isolation of Group B streptococci

Package contents:

Item code (packaging) REF	Description	Primary packaging code:	Number of products
PRM1840V20	Substrate poured into petri dishes of ø90	PRM1840	20
PRM1840V60			60
PRM1840V240			240
PRM1840M40	Substrate poured into petri dishes of ø50		40

Directions

Surface spread the test inoculum aseptically on the plate. After incubation, observe growth and color of colonies.

Principle and interpretation

Group B *Streptococcus* is a leading infection causing illness and death in newborns. Group B streptococci can also cause serious diseases in pregnant women, the elderly, and adults with other illnesses. GBS normally reside in the vagina of women and rectum of men and women (1). In newborns, group B strep is the most common cause of sepsis (infection of the bloodstream) and meningitis (infection of the lining and fluid surrounding the brain) and a common cause of pneumonia.

In adults, group B strep can rarely lead to serious bloodstream infections, urinary tract infections, skin infections, and pneumonia, especially in people with weak immune systems. Heavy colonization of the maternal genital tract is associated with colonization of infants and risk of neonatal disease (2).

The sample collection is usually done by collection of vaginal and rectal swab between 35 and 37 weeks of pregnancy. The swab is then processed on Chromogenic Strep B Selective Agar Plate. For the conventional methods optimum recovery is however achieved by selective enrichment into Todd Hewitt broth with colistin and nalidixic acid and then subculture on Blood Agar (3,4).

Protein hydrolysate provides essential nutrients for the growth of Streptococci. Buffers present provides buffering to the medium. Selective agents in the medium inhibits accompanying flora. One of the substrate in the chromogenic mixture is cleaved by beta glucosidase possessed by Group B Streptococci resulting in blue coloured colonies.

Quality control

The data and results of quality control are given in the certificate of analysis for each lot.

Storage and shelf life

Storage between 2-8°C. Use before expiry date on the label.


Warning and precautions

In vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Symbols used on labels

	European Conformity mark		This side up
	is an in vitro diagnostic medical device (IVD)		Catalogue number
	Do not expose directly to sunlight		Batch code
	Consult instructions for use		Use-by date
	Do not re-use		Temperature limit
	Pack size		Manufacturer
	European Authorized Representative (Authorised Representative)		

	Salus Cons kft. 6722 Szeged, Bécsi krt 23, HUNGARY e-mail: office@saluscons.com
---	--

Reference

1. Anthony BF, Okada DM, Hobel CJ. Epidemiology of group B Streptococcus: longitudinal observations during pregnancy. J Infect Dis 1978; 137:524-30.
2. Murray P.R., Baron J.H., Manual of Clinical Microbiology Murray P. R., Baron J. H., Tenover F.C., Tenover F.C., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
3. Prevention of perinatal group B Streptococcal disease: a public health perspective. Centres for Disease control and Prevention. MMWR Recomm Rep 1996; 51:1-22
4. NHS Processing swabs for Group B Streptococcal carriage Issue no.2.1,2006