UPUTSTVO ZA UPOTREBU (SRB)

Thayer Martin Agar Plate w/VCNT

Podloga za selektivnu izolaciju gonococci iz patoloških uzoraka.

Sadržaj pakovanja:

Šifra artikla (pakovanja) REF	Opis	Šifra primarnog pakovanja:	Broj podloga
PRM413V20	Podloga izlivena u petri posudama od ø90	PRM413	20
PRM413V60			60
PRM413V240			240
PRM413M40	Podloga izlivena u petri posudama od ø50		40

Uputstva

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
INSTRUCTION FOR USE (EN)

Thayer Martin Agar Plate w/VCNT

Medium is used for selective isolation of Gonococci from pathological specimens.

Package contents:

Item code (packaging) REF	Description	Primary packaging code:	Number of products
PRM413V20	Substrate poured into petri dishes of ø90	PRM413	20
PRM413V60			60
PRM413V240			240
PRM413M40	Substrate poured into petri dishes of ø50		40

Directions

Surface spread the test inoculum aseptically on the plate.

Principle and interpretation

Carpenter and Morton reported an improved medium to isolate Gonococci in 24 hours (1). Later on the efficiency of GC medium supplemented with haemoglobin and yeast concentrate was demonstrated for isolating gonococci (2). Subsequently Thayer and Martin Medium was developed for the primary isolation of *Neisseria gonorrhoeae* and *Neisseria meningitidis* from specimens containing mixed flora collected from throat, vagina, rectum and urethra (3, 4). Thayer and Martin (4) used Vancomycin, Colistin and Nystatin. Martin and Lester (5) used an additional antibiotic Trimethoprim to make the medium selective. Special peptone provides nutrients to the organisms while starch neutralizes the toxic fatty acids if present in the agar. Haemoglobin provides the X factor whereas the V factor (NAD) is provided by the added supplement. Medium also supplies vitamins, amino acids, coenzymes etc. which enhances the growth of pathogenic *Neisseria* spp. Vancomycin and colistin inhibits Gram-positive and Gram-negative bacteria respectively (6). Nystatin inhibits fungi. This medium may inhibit *Haemophilus* species. Some strains of *Capnocytophaga* species may grow on this medium when inoculated with oropharyngeal specimens.

Quality control

The data and results of quality control are given in the certificate of analysis for each lot.

Storage and shelf life

Storage between 2-8°C. Use before expiry date on the label.



Warning and precautions


In vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Symbols used on labels

	European Conformity mark		This side up
	is an in vitro diagnostic medical device (IVD)		Catalogue number
	Do not expose directly to sunlight		Batch code
	Consult instructions for use		Use-by date
	Do not re-use		Temperature limit
	Pack size		Manufacturer
	European Authorized Representative (Authorised Representative)		

	Salus Cons kft. 6722 Szeged, Bécsi krt 23, HUNGARY e-mail: office@saluscons.com
---	--

Reference

1. Carpenter and Morton, 1947, Proc. N.Y. State Assoc. Public Hlth. Labs., 27:58.
2. Carpenter et al, 1949, Am. J. Syphil. Gonorrh. Vener. Dis., 33:164.
3. Martin, Billings, Hackney and Thayer, 1967, Public Hlth. Rep., 82:361.
4. Thayer J. and Martin J.E. Jr., 1966, Public Health Rep., 81:559.
5. Martin J.E. Jr. and Lester A., 1971, HSMHA Hlth. Service Rep., 86(1):30.
6. MacFaddin J., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.