UPUTSTVO ZA UPOTREBU (SRB)

Sabouraud Dextrose Agar Plate w/Chloramphenicol (50mg/L) and Cycloheximide (500mg/L)

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
INSTRUCTION FOR USE (EN)

Sabouraud Dextrose Agar Plate w/Chloramphenicol and Cycloheximide

Medium is recommended for selective cultivation of yeasts and moulds.

Package contents:

Item code (packaging) REF	Description	Primary packaging code:	Number of products
PRM664V20	Substrate poured into petri dishes of ø90	PRM664	20
PRM664V60			60
PRM664V240			240
PRM664M40	Substrate poured into petri dishes of ø50		40

Directions

Surface spread the test inoculum aseptically on the plate. After incubation, observe growth and color of colonies.

Principle and interpretation

Sabouraud Dextrose Agar was originally formulated by Sabouraud (1) and further modified by Emmons (2) by reducing dextrose content and adjusting the pH close to neutral.

Peptic digest of animal tissue is the source of nitrogenous growth factors while dextrose provides an energy source for the growth of microorganisms. The media can be rendered selective for fungi by antibiotics such as Chloramphenicol (4) and Cycloheximide (5), which inhibit some bacteria as well as some saprophytic and pathogenic fungi. This medium inhibits fungi like *Cryptococcus neoformans*, *Aspergillus*, *Nocardia*, certain *Candida* species but allow the dermatophytes to grow well.

Quality control

The data and results of quality control are given in the certificate of analysis for each lot.

Storage and shelf life

Storage between 2-8°C. Use before expiry date on the label.


Warning and precautions


In vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Symbols used on labels

	European Conformity mark		This side up
	is an in vitro diagnostic medical device (IVD)		Catalogue number
	Do not expose directly to sunlight		Batch code
	Consult instructions for use		Use-by date
	Do not re-use		Temperature limit
	Pack size		Manufacturer
	European Authorized Representative (Authorised Representative)		

	Salus Cons kft. 6722 Szeged, Bécsi krt 23, HUNGARY e-mail: office@saluscons.com
---	--

Reference

- Sabouraud R., 1892, Ann. Dermatol. Syphilol., 3:1061.
- Emmons C., Binford C., Uty J. and Kwon-Chung, 1970, Medical Mycology, 2nd ed., Philadelphia: Lea and Febiger.
- Diagnostic Procedures, 1963, 4th ed., APHA
- Ajello L., 1957, J. Chron. Dis., 5:545.,
- MacFaddin J. F., 1985, Media For Isolation-Cultivation Identification - Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore.